Table of contents

1 alarm	Page number: 2
2 bedsidemonitor	Page number: 3
3 bloodpressure	Page number: 4
4 breathingrate	Page number: 5
5 centralstation	Page number: 6
6 nurse	Page number: 7
7 patient	Page number: 8
8 pulserate	Page number: 9
9 reading	Page number: 10
10 role	Page number: 11
11 temperature	Page number: 12
12 user_account	Page number: 13
13 Relational schema	Page number: 14

Page number: 1/14

1 alarm

Creation: Nov 20, 2019 at 01:27 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
reading_id	int(11)		No			-> reading.id ON UPDATE CASCADE ON DELETE CASCADE		
alarmLevel	varchar(20)		No					
dateTimeTri gger	datetime		No					
dateTimeMu ted	datetime		No					

2 bedsidemonitor

Creation: Nov 20, 2019 at 01:25 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
patient_id	int(11)		No			-> patient.id ON UPDATE CASCADE ON DELETE CASCADE		
centralStati on_id	int(11)		No			-> centralstation.id ON UPDATE CASCADE ON DELETE CASCADE		

3 bloodpressure

Creation: Nov 20, 2019 at 01:21 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
patient_id	int(11)		No			-> patient.id ON UPDATE CASCADE ON DELETE CASCADE		
bloodPressu reValue	double		No					
bloodPressu reTime	time		No					
bloodPressu reDate	date		No					

4 breathingrate

Creation: Nov 20, 2019 at 12:40 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
patient_id	int(11)		No			-> patient.id ON UPDATE CASCADE ON DELETE CASCADE		
breathingRa teValue	double		No					
breathingRa teTime	time		No					
breathingRa teDate	date		No					

5 centralstation

Creation: Nov 20, 2019 at 12:44 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
alarm_id	int(11)		No			-> alarm.id ON UPDATE CASCADE ON DELETE CASCADE		

6 nurse

Creation: Nov 20, 2019 at 01:13 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
nurse_id	int(11)		No					

Page number: 7/14

7 patient

Creation: Nov 20, 2019 at 12:47 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
reading_id	int(11)		No			-> reading.id ON UPDATE CASCADE ON DELETE CASCADE		
fullName	varchar(20)		No					
NRIC	int(30)		No					
address	varchar(100		No					
email	varchar(20)		No					
gender	varchar(10)		No					
phone	int(30)		No					
age	int(11)		No					

8 pulserate

Creation: Nov 20, 2019 at 12:49 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
patient_id	int(11)		No			-> patient.id ON UPDATE CASCADE ON DELETE CASCADE		
pulseRateVa lue	double		No					
pulseRateTi me	time		No					
pulseRateD ate	date		No					

9 reading

Creation: Nov 20, 2019 at 12:50 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
minValueBlo odPressure	double		No					
maxValueBl oodPressure	double		No					
minValueTe mperature	double		No					
maxValueTe mperature	double		No					
minValueBr eathingRate	double		No					
maxValueBr eathingRate	double		No					
minValuePul seRate	double		No					
maxValuePu IseRate	double		No					

10 role

Creation: Nov 20, 2019 at 12:52 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
position	varchar(20)		No					
user_accoun t_id	int(11)		No			-> user_account.id ON UPDATE CASCADE ON DELETE CASCADE		

Page number: 11/14

11 temperature

Creation: Nov 20, 2019 at 12:53 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
patient_id	int(11)		No			-> patient.id ON UPDATE CASCADE ON DELETE CASCADE		
temperature Value	double		No					
temperature Time	time		No					
temperature Date	date		No					

Page number: 12/14

12 user account

Creation: Nov 20, 2019 at 12:55 PM

Column	Туре	Attributes	Null	Default	Extra	Links to	Comments	MIME
id	int(11)		No		auto_inc rement			
password	varchar(20)		No					
loginStatus	varchar(20)		No					
NRIC	int(30)		No					
address	varchar(100)		No					
phone	int(20)		No					
dateTimeRe gister	datetime		No					
dateTimeLo gin	datetime		No					
dateTimeLo gOut	datetime		No					

Page number: 13/14

